

PowerNurse™

TRANSFER OPTIONS

PowerNurse can transfer patients utilizing three bedding options:

1 A chuck under the patient's midsection and pillows placed under the knees and head.



2 A temporary transfer pad bed-rolled under the patient.



3 A custom two-part mattress designed for indefinite use under the patient.



- Options 1 and 3 do NOT require a bed roll and therefore can be accomplished by a SINGLE NURSE in a ZERO LIFT environment.
- Option 2 provides as much patient comfort as option 3, but involves bed rolling the patient.
- All three transfer options can be used, depending on the need or circumstance.



Hinged transfer pad allows full bed coverage and fits narrow stretchers.

The solution to your lateral transfer challenges

PowerNurse™

SPECIFICATIONS

SIZE: 77" x 28" x 2.25"
WEIGHT: < 100 lbs.
SPEED: 1 & 2 in/sec.
POWER: 110-240V, 50-60 Hz




1679 Lowell Road
Concord, MA 01742
508.951.8842
www.astirtechnologies.com



PowerNurse™
Single Nurse Operated
Powered Patient Transfer Device

Astir Technologies introduces *PowerNurse*, the most important breakthrough in patient transfer since the slide board.



*PowerNurse*TM

ADVANTAGES

- Zero lift device*
- Single nurse operated*
- 3-minute transfer*
- Increase O. R. throughput*
- Decrease nurse MSD injuries*
- Reduce staff turnover*
- Reduce overtime and replacement staff*
- Decrease hospital costs*



*PowerNurse*TM

The solution to your lateral transfer challenges



www.astirtechnologies.com

BRAKES AND RAILING "LOCK-OUT"—
reduces errors and enhances ease of use

MOVE COMMAND—
hold down to move or rotate in pre-selected mode

BELT CONTROL—
select from 3 modes

- burrow
- pad alignment
- transfer

